




- ▶ Mobile disinfection column including disinfection dispenser and drip tray
- ✓ Column made from stainless steel, height 1.17 m
- ✓ Dispenser with spray system, content 0.9 litres
- ✓ Drip tray with removable collection tray




- ▶ Stable on flat surfaces thanks to robust stainless steel design



- ▶ Flexible use
- ✓ No attachment to wall or floor required



- ▶ Clean as well as hygienic
- ✓ Thanks to the drip tray, the floor stays dry



Disinfection dispenser, drip tray, and disinfection column as a mobile hygiene package. Thanks to the robust stainless steel design, the set can be used as required, without fixation, anywhere on an even surface.

- | | | | |
|------------------------------|---|-----------|---------|
| • Set comprises: | Disinfection column D1 1170 (construction kit: 1 pedestal, 1 column)
Disinfectant dispenser PS 0.9L-W
Drip tray PS 0.9L-W | • Weight: | 8.53 kg |
| • Content: | 0,9 litres | | |
| • Cartridge-free: | Yes | | |
| • Can be refilled: | Yes | | |
| • Operation: | One-hand operation
Elbow operation
Manual | | |
| • Filling level check: | Transparent plastic container | | |
| • Collection tray: | Detachable
Dishwasher-suitable
Size: W 147 x D 240 x H 40 mm | | |
| • Delivery state: | Building kit (simple assembly) | | |
| • Properties: | - | | |
| • Including: | Attachment material
1 wall holder | | |
| • Not included in delivery: | Filling | | |
| • Material: | Stainless steel
Aluminium
Polypropylene
Plastic ABS | | |
| • Suitable for Euro bottles: | 1000 ml | | |
| • Dosing: | Pump system
Spray system
Up to 3 ml per stroke | | |
| • Important information: | - | | |
| • Size: | W 300 x D 400 x H 1,170 mm | | |

SKINTASTIC LEOCID® SEPT P7



- Order quantity unit: 1 box (12 bottles)
- Content: 12 x 1 litre
- Model: Liquid
- Ambient temperature: 5 °C to 25 °C
- Ph level: 7
- Properties: Partially viricidal
Bactericidal, yeasticidal, tuberculocidal
Perfume-free
- Designed for: Disinfectant dispenser
- Application: Allow 3 ml of undiluted liquid to act for 30 seconds (as recommended by RKI)

Code-No. 173084
GTIN 4029888031308